

PMW CONTROL SWAB

Device Name

Product Name: PMW-SWAB

Model Name: PMW CONTROL SWAB | PMW CONTROL SWAB L15 | PMW CONTROL SWAB-B

General Description

The PMW CONTROL SWAB, PMW CONTROL SWABL15 and PMW CONTROL SWAB-B are feature POLYURETHANE fiber swabs that optimize specimen collection. The swabs also feature a molded breakpoint that allows you to safely and easily break off the swab stick, and several breakpoint options are available.

Intended Purpose

A nasal or oral swab is a method for collecting a clinical test sample of nasal cavity or oral cavity secretions from the nose and throat. The sample is then analyzed for the presence of organisms or other clinical markers for disease.

Performance Characteristics

PMW-SWAB (PMW CONTROL SWAB, PMW CONTROL SWAB L15 and PMW CONTROL SWAB-B) is a specimen collection swab that has polyurethane fiber. The polyurethane specimen collection swabs have a solid plastic shaft with a molded breakpoint site.

To collect the sample, the swab is inserted in the nasal and oral cavity. The swab is then rotated for a specified period time to collect secretions, then the swab is removed.

PMW-SWAB is a method for collecting a clinical test sample of nasal or oral secretions. The sample is then analyzed for the presence of organisms or other markers for disease.

Indications

The swab is suitable for short contact with nasal and oral cavity of the human body.

Contra-indications

There are no known specified contraindications.

Patient Target Group

Age	Not relevant	Weight	Not relevant	Nationality	Not relevant
Gender	Not relevant	Health	Not relevant	Patient State	Not relevant

Intended Users

This medical device is intended to be used only by healthcare professionals such as doctors and medical personnel (e.g., research nurse or physician)

Residual Risks

There is no residual risk with this medical device after conducting risk management activities.

Precautions and warnings

- This product is disposable and reuse is strictly prohibited.
- This product should not be used if the package is damaged or opened.
- Do not use if expired expiration date.
- It should be used only for intended purpose.
- Do not use the device perpendicularly. Please use the device in an oblique direction.

PMW CONTROL SWAB

Attention

- Take care as this may break if excessive force is applied during specimen collection.
- If this product is exposed to chemical and physical sterilization or micro-biostatic process, its intended performance may not be achieved.
- User should fully aware of usage before use.
- When handling specimen, be sure to have protective equipment to avoid any risk of infection and tested in accordance with bio-safety guidelines.
- Do not use unless peel-pouch is intact.
- Please be aware that the head part should not be touched.
- It should be disposed in a disposal box after used.
- Please be aware of using the device perpendicularly, the break point of the device may break during use if the user uses the device perpendicularly.

Undesirable Side-effects

The tip of the device may break during use, and the patient may swallow the broken part with the specimen. While collecting the specimen by using a swab, it can break and stab the nose and throat.

Preparatory Treatment or Handling

- Storage at room temperature(1~30°C) as of sealed status.
- Check the sealed status if damaged or abnormalities and product damage before use and do not use if the package is damaged or opened.
- Check the expiration date of the product before use and do not use if expired.
- User should fully understand how to use this product.
- Wear gloves to prevent any infection.

Shelf-life (Life time)

- 3 years from the date of manufacture

Single Use

- This medical device is intended for single use.
- Information on known characteristics and technical factors known to the manufacturer that could pose a risk if the device were to be re-used
- There is a risk of infection if the device is reused.

Information for Safe Disposal

- This device is disposable for single use only.
- Once used, it must be disposed of in accordance with the regulations, and reuse is prohibited.
- Wear gloves to prevent any infection.
- It should be disposed in a disposal box after used.

Notice to Users/Patients


Any serious incident that has occurred in relation to this medical device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Product Components

No.	Model	Contents
1	PMW CONTROL SWAB	A swab that is intended to collect specimens from the nasal and/or oral is included in one blister packaging
2	PMW CONTROL SWAB-B	
3	PMW CONTROL SWAB L15	

PMW CONTROL SWAB

How to Use



- ① Peel off the sealed film package then, pull the product out after grab the body handle.
- ② Collecting specimen from the mouth, nasal cavity, ear or anus. In order to reduce the possibility of other contamination, the TIP part of the head is in contact with the suspected infection area quickly.
- ③ Rotate the body handle a few times to ensure that the specimen are effectively collected to the head end.

Symbols of Packaging

	CE marking		Medical Device		Manufacturer
	Authorized representative in the European Community		Date of manufacturer		Use-by date
	Lot Number		Catalogue Number		Do not use if package is damaged
	Temperature limit		Do not re-use		Caution



PMW corp.

Headquarters : Room 904, 48, Centum jungang-ro, Haeundae-gu, Busan, Republic of Korea

+82-51-507-6364 / info@pmwcorp.co.kr / info@pmwglobal.co.kr

Site Address : Room A, 2F, 135, Cheoyongsaneop-3gil, Cheongnyang-eup, Ulju-gun, Ulsan, Republic of Korea

+82-52-227-6364 / info@pmwcorp.co.kr / info@pmwglobal.co.kr



CMC Medical Devices & Drugs S.L.

C/Horacio Lengo N° 18, CP 29006, Mlaga, Spain / +34-951-214-054 / info@cmcmedicaldevices.com